


power one ACCU *plus*.

Rechargeable hearing aid batteries.



A quality feature of **power one ACCU *plus*** hearing aid batteries is the unique nano surface coating of the stainless steel battery casing.

The **power one p675accu** is characterized by the following features

- Nano surface coating to ensure great corrosion resistance, even if exposed to large quantities of sweat
- High stability and unlimited function in a wide range of ambient conditions
- Environmentally-friendly technology
- 0% mercury, lead, and cadmium
- Made in Germany.

Please note that original **power one** batteries can be identified by their coloured edging.

Call us on +49 (0)79 61 921 790 or
contact your distribution partner.

www.powerone-batteries.com

power one 675charger.



The **power one 675charger** is for charging one or two rechargeable **power one ACCU plus** type p675accu hearing aid batteries.

Flat **power one** p675accu batteries are fully charged within five hours. To improve their safety and durability, batteries that are not completely flat are automatically discharged before the charging process.

Please only use **power one ACCU plus** rechargeable batteries in the **power one 675charger**.



Shown original size

www.powerone-batteries.com